

Anti-POMC antibody (1-267 aa) (STJ11106269)

STJ11106269

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	5443
Gene Symbol	POMC
Uniprot ID	COLI_HUMAN
Immunogen Region	1-267 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-267 of human POMC (NP_000930.1).
Immunogen Sequence	MPRSCCSRSGALLALLLQA SMEVRGWCLESSQCQDLTTE SNLLECIRACKPDLSAETPM FPGNGDEQPLTENPRKYVMG HFRWDRFGRNRNSSSGSSGA GQKREDVSAGEDCGPLPEGG PEPRSDGAKPGPREGKRYSYS MEHFRWGKPVGKKRRPVKVY PNGAEDESAAEFLEFKREL TGQRLREGDGPDPADDGAG AQADLEHSLLVAAEKKDEGP YRMEHFRWGSPPKDKRYGG