

Anti-APBA2 antibody (400-500 aa) (STJ11106271)

STJ11106271

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	321
Gene Symbol	APBA2
Uniprot ID	APBA2_HUMAN
Immunogen	
Immunogen	400-500 aa
Region	
Specificity	A synthetic peptide corresponding to a sequence within amino acids 400-500 of human APBA2 (NP_005494.2).
Immunogen	AVSRVKRMQKAAKIKKKANS EGDAQTLTEVDLFISTQRIK VLNADTQETMMDHALRTISY IADIGNIVVLMARRRMPRSA
Sequence	SQDCIETTPGAQEGKKQYKM I