

Anti-RAMP2 antibody (STJ11106283)

STJ11106283

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	10266
Gene Symbol	RAMP2
Uniprot ID	RAMP2_HUMAN
Immunogen	
Immunogen Region	
Specificity	A synthetic peptide of human RAMP2.
Immunogen Sequence	