

Anti-NID2 antibody (900-1200 aa) (STJ11106287)

STJ11106287

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	22795
Gene Symbol	NID2
Uniprot ID	NID2_HUMAN
Immunogen	
Immunogen Region	900-1200 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 900-1200 of human NID2 (NP_031387.3).
Immunogen Sequence	RCHPAATCYNTPGSFSCRCQ PGYYGDGFQCIPDSTSSLTP CEQQQRHAQAQYAYPGARFH IPQCDEQGNFLPLQCHGSTG FCWCVDPDGHEVPGTQTTPPG STPPHCGPSPEPTQRPPTIC ERWRENLEHYGGTPRDDQY VPQCDDLGHFIPLQCHGKSD FCWCVDKDGREVQGTRSQPG TTPACIPTVAPPMVRPTPRP DVTPPSVGTFLLYTQGQIG YLPLNGTRLQKDAAKTLLS