

Anti-TRAPPC8 antibody (150-450 aa) (STJ11106288)

STJ11106288

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	22878
Gene Symbol	TRAPPC8
Uniprot ID	TPPC8_HUMAN
Immunogen	
Immunogen Region	150-450 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 150-450 of human TRAPPC8 (NP_055754.2).
Immunogen Sequence	LVASSSEAEPVQFSKLSQE QHRIQHNSDYSYPKWFIPNT LKYYVLLHDVSAGDEQRAES IYEEMKQKYGTQGCYLLKIN SRTSNRASDEQIPDPWSQYL QKNSIQNQESYEDGPGCTITS NKNSDNLLSLDGLDNEVKD GLPNNFRAHPLQLEQSSDPS NSIDGPDHLRSASSLHETKK GNTGIIHGACLTLDHDIRI QFIQFTFRGLLPHEKTIR QLNDQLISRKGLSRSLFSA