

Anti-CBX7 antibody (1-110 aa) (STJ11106294)

STJ11106294

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	23492
Gene Symbol	CBX7
Uniprot ID	CBX7_HUMAN
Immunogen	
Immunogen Region	1-110 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human CBX7 (NP_783640.1).
Immunogen Sequence	MELSAIGEQVFAVESIRKKR VRKGKVEYLVKWKGWPPKYS TWEPEEHILDPRLMAYEEK EERDRASGYRKRGPKPKRLL LQRLYSMDLRSSHKAGKEK LCFSLTCPLG