

**Anti-FOXP3 antibody (1-70 aa) (STJ11106301)**

STJ11106301

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	WB/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	50943
<b>Gene Symbol</b>	FOXP3
<b>Uniprot ID</b>	FOXP3_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	1-70 aa
<b>Specificity</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human FOXP3 (NP_054728.2).
<b>Immunogen Sequence</b>	MPNPRPGKPSAPSLALGPSP GASPSWRAAPKASDLLGARG PGGTFQGRDLRGGAHASSSS LNPMPPSQLQ

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081