

Anti-RAPGEF6 antibody (145-240 aa) (STJ11106303)

STJ11106303

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	51735
Gene Symbol	RAPGEF6
Uniprot ID	RPGF6_HUMAN
Immunogen	
Immunogen Region	145-240 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 145-240 of human RAPGEF6 (NP_001157858).
Immunogen Sequence	INYGGERQTITDDDEVNSYL SLPADLTkMHLTENPHQVt HVSSSQSGCSIASDSGSSSL SDIYQATESEVGdVDLTRLp EGPVdSEdEEEEdEEI