

Anti-INO80 antibody (1400 aa) (STJ11106304)

STJ11106304

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	54617
Gene Symbol	INO80
Uniprot ID	INO80_HUMAN
Immunogen	
Immunogen Region	1400 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1400 to the C-terminus of human INO80 (NP_060023.1).
Immunogen Sequence	PASITGSVSDTVNGISIQEM PAAGRGHSAERSRGRPKGSGS TAKGAGKGRSRKSTAGSAAA MAGAKAGAAAAAAYG YNVSKGISASSPLQTSVLRP AGLADFGPSSASSPLSSPLS KGNNVPGNPKNLHMTSSLAP DSLVRKQGKGTNPSSGGR