

Anti-FAIM antibody (1-179 aa) (STJ11106308)

STJ11106308

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	55179
Gene Symbol	FAIM
Uniprot ID	FAIM1_HUMAN
Immunogen Region	1-179 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-179 of human FAIM (NP_060617.1).
Immunogen Sequence	MTDLVAVWDVALSDGVHKIE FEHGTTSQKRVVYVDGKEEI RKEWMFKLVGKETFYVGA AK TKATINIDAISGFAYEYTL E INGKSLKKYMEDRSKTTNTW VLHMDGENFRIVLEKDAMDV WCNGKKLETAGEFVDDGTET HFSIGNHDCYIKAVSSGKRK EGIIHTLIVDNREIPEIAS