

Anti-SYT12 antibody (95-200 aa) (STJ11106318)

STJ11106318

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	91683
Gene Symbol	SYT12
Uniprot ID	SYT12_HUMAN
Immunogen	
Immunogen Region	95-200 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 95-200 of human SYT12 (NP_001171351.1).
Immunogen Sequence	KGSLSIEDTFESISELGP LMGRELDAPYGTLRKSQSA DSLNSISSVSNTFGQDFTLG QVEVSMEYDTASHTLNVAVM QGKDLLEREASFESCFMRV SLLPDE