

Anti-CADPS2 antibody (1-80 aa) (STJ11106319)

STJ11106319

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	93664
Gene Symbol	CADPS2
Uniprot ID	CAPS2_HUMAN
Immunogen	
Immunogen Region	1-80 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human CADPS2 (NP_001009571.2).
Immunogen Sequence	MLDPSSEEESEGLEEEESR DVLVAAGSSQRAPPAPTREG RRDAPGRAGGGGAARSVSPS PSQLSEGRDEPQRQLDDEQE