

Anti-EXOSC6 antibody (150-250 aa) (STJ11106320)

STJ11106320

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	118460
Gene Symbol	EXOSC6
Uniprot ID	EXOS6_HUMAN
Immunogen	
Immunogen Region	150-250 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human EXOSC6 (NP_478126.1).
Immunogen Sequence	LLEDGGSALAAALTAALAL ADAGVEMYDLVGCGLSLAP GPAPTWLLDPTRLEEERAAA GLTVALMPVLNQVAGLLGSG EGGLTESWAEAVRLGLEGCQ R