

Anti-IRX2 antibody (200-300 aa) (STJ11106323)

STJ11106323

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	153572
Gene Symbol	IRX2
Uniprot ID	IRX2_HUMAN
Immunogen	
Immunogen Region	200-300 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human IRX2 (NP_150366.1).
Immunogen Sequence	KDESPDKAQEGTETSAEDEG ISLHVDSLTDHSCSAESDGE KLPCRAGDPLCESGSECKDK YDDLEDDDDDEEGERGLAP PKPVTSSPLTGLEAPLLSPP P