

Anti-SLC35F1 antibody (1-100 aa) (STJ11106327)

STJ11106327

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	222553
Gene Symbol	SLC35F1
Uniprot ID	S35F1_HUMAN
Immunogen	
Immunogen Region	1-100 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SLC35F1 (NP_001025029.2).
Immunogen Sequence	MIPPEQPQQQLQPPSPAPPN HVVTTIENLPAEGSGGGGSL SASSRAGVRQRIKVLNREM LISVALGQVLSLLICIGILT SKYLSEDFHANTPVFQSFLN