

Anti-MIOX antibody (24-431 aa) (STJ11106330)

STJ11106330

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	431705
Gene Symbol	ASTL
Uniprot ID	ASTL_HUMAN
Immunogen	
Immunogen Region	24-431 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 24-431 of human MIOX (NP_001002036.3).
Immunogen Sequence	APLASSCAGACGTSFPDGLT PEGTQASGDKDIPAINQGLI LEETPESSFLIEGDIIRPSP FRLLSATSNKWPMGGSGWVE VPFLSSKYDEPSRQVILEA LAEFERSTCIRFVTYQDQRD FISIIPMYGCFSSVGRSGGM QVSLAPTCLQKGRGIVLHE LMHVLGFWHEHTRADRDRYI RVNWNEILPGFEINFIKSQS SNMLTPYDYSSVMHYGRLAF SRRGLPTITPLWAPSVHIG