

Anti-SETD1A antibody (300-400 aa) (STJ11106335)

STJ11106335

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	9739
Gene Symbol	SETD1A
Uniprot ID	SET1A_HUMAN
Immunogen	
Immunogen Region	300-400 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human SETD1A (NP_055527.1).
Immunogen Sequence	RSNSYQDAFSRRHFSASSA STTASTAIAATTAATASSA SSSSLSSSSSSSSSSSQF RSSDANYPAYYESWNRVQRH TSYPPRRATREEPPGAPFAE N