

Anti-SIDT2 antibody (19-290 aa) (STJ11106336)
STJ11106336

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	51092
Gene Symbol	SIDT2
Uniprot ID	SIDT2_HUMAN
Immunogen Region	19-290 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 19-290 of human SIDT2 (NP_001035545.1).
Immunogen Sequence	HLGVLGPKNVSQKDAEFERT YVDEVNSELVNIYTFNHTVT RNRTEGVRVSVNVLNKQKGA PLLFVVRQKEAVVSFQVPLI LRGMFQRKYLYQKVERTLCQ PPTKNESEIQFFYVDVSTLS PVNTTYQLRVSRMDDFVLRT GEQFSFNTTAAQPQYFKYEF PEGVDSVIVKVTSNKAFPCS VISIQDVLCPVYDLNNAVAF IGMYTMTKKAITVQRKDF PSNSFYVVVVKTEDQACG