

Anti-GBP5 antibody (477-586 aa) (STJ11106341)

STJ11106341

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	115362
Gene Symbol	GBP5
Uniprot ID	GBP5_HUMAN
Immunogen	
Immunogen Region	477-586 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 477-586 of human GBP5 (NP_443174.1).
Immunogen Sequence	DQALTETEKKKKEAQVKAEA EKAEAQRLAAIQRQNEQMMQ ERERLHQEQVRQMEIAKQNW LAEQQKMQEQQMQEQAAQLS TTFQAQNRSLSELQHAQRT VNNDPCVLL