

Anti-Acetyl-FOXO1-K262/K265/K274 antibody (221-320 aa) (STJ11106349)

STJ11106349

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:100-1:500
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	2308
Gene Symbol	FOXO1
Uniprot ID	FOXO1_HUMAN
Immunogen	
Immunogen Region	221-320 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 221-320 of human FOXO1 (NP_002006.2).
Immunogen Sequence	SKFIRVQNEGTGKSSWWMLN PEGGKSGKSPRRRAASMDNN SKFAKSRRAAKKASLQSG QEGAGDSPGSQFSKWPA SPG SHSNDDFDNWSTFRPRTSSN ASTISGRLSPIMTEQDDLGE GDVHSMVYPPSAAKMASTLP SLSEISNPENMENLLDNLNL LSSPTSLTVSTQSSPGTMMQ QTPCYSFAPPNTSLNSPSPN YQKYTYGQSSMSPLPQMPIQ TLQDNKSSYGMSQYNCAP



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081