

Anti-COL12A1 antibody (24-270 aa) (STJ11106353)

STJ11106353

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	1303
Gene Symbol	COL12A1
Uniprot ID	COCA1_HUMAN
Immunogen Region	24-270 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 24-270 of human COL12A1 (NP_004361.3).
Immunogen Sequence	EVDPPSDLNFKIIDENTVHM SWAKPVDPIVGYRITVDPTT DGPTEFTLSASTTETLLSE LVPETEVVVTITSYDEVEES VPVIGQLTIQTGSSTKPVEK KPGKTEIQKCSVSAWTDLVF LVDGSWSVGRNNFKYILDFI AALVSAFDIGEEKTRVGWVQ YSSDTRTEFNLNQYYQRDEL LAAIKKIPYKGGNTMTGDAI DYLVKNFTTESAGARVGFPK VAIITDGKSQDEVEIPAR