

Anti-DAG1 antibody (776-895 aa) (STJ11106355)

STJ11106355

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	1605
Gene Symbol	DAG1
Uniprot ID	DAG1_HUMAN
Immunogen	
Immunogen Region	776-895 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 776-895 of human DAG1 (NP_004384.4).
Immunogen Sequence	RKKRKGKLTLEDQATFIKKG VPIIFADELDDSKPPPSM PLILQEEKAPLPPEYPNQS VPETPLNQDTMG EYTPLRD EDPNAPPYQPPPPFTAPMEG KGSRPKNMTPYRSPPPYVPP