

Anti-GK2 antibody (1-150 aa) (STJ11106359)

STJ11106359

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	2712
Gene Symbol	GK2
Uniprot ID	GLPK2_HUMAN
Immunogen	
Immunogen Region	1-150 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human GK2 (NP_149991.2).
Immunogen Sequence	MAAPKTAAVGPLVGAVVQGT NSTRFLVFSKTAELLSHHK VELTQEFKPEGWVEQDPKEI LQSVYECIARTCEKLDLNI DISNIKAVGVSNQRETTVIW DKLTGEPLYNAVWLDLRTQ TTVEDLSKKIPGNSNFVKSK TGLPLSTYFS