

Anti-HMGN1 antibody (1-100 aa) (STJ11106363)

STJ11106363

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	3150
Gene Symbol	HMGN1
Uniprot ID	HMGN1_HUMAN
Immunogen	
Immunogen Region	1-100 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HMGN1 (NP_004956.5).
Immunogen Sequence	MPKRKVSSAEGAAKEEPKRR SARLSAKPPAKVEAKPKKAA AKDKSSDKKVQTKGKRGAKG KQAEVANQETKEDLPAENGE TKTEESPASDEAGEKEAKSD