

Anti-HTR5A antibody (220-285 aa) (STJ11106364)

STJ11106364

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	3361
Gene Symbol	HTR5A
Uniprot ID	5HT5A_HUMAN
Immunogen	
Immunogen Region	220-285 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 220-285 of human HTR5A (NP_076917.1).
Immunogen Sequence	YWKIYKAAKFRVGSRKTNV SPISEAVEVKDSAKQPMVF TVRHATVTFQPEGDTWREQK EQRAAL