

Anti-ZBTB14 antibody (1-280 aa) (STJ11106368)

STJ11106368

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	7541
Gene Symbol	ZBTB14
Uniprot ID	ZBT14_HUMAN
Immunogen	
Immunogen Region	1-280 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ZBTB14 (NP_003400.2).
Immunogen Sequence	MEFFISMSETIKYNDHKT LFLKTLNEQRLEGEFCDAI VEDVKFRAHRCVLAACSTY FKFLFKKLEVDSSSVIEIDF LRSDIFEEVLNYMYTAKISV KKEDVNLMMSSGQILGIRFL DKLCSQKRDVSSPDENNGQS KSKYCLKINRPIGDAADTQD DDVEIGDQDDSPSDDTVEG TPPSQEDGKSPTTTTLRVQEA ILKELGSEEVKVNVCYQGEV ESMETPESKDLGSQTPQAL