

Anti-MICU1 antibody (50-200 aa) (STJ11106379)

STJ11106379

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	10367
Gene Symbol	MICU1
Uniprot ID	MICU1_HUMAN
Immunogen	
Immunogen Region	50-200 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 50-200 of human MICU1 (NP_001182447.1).
Immunogen Sequence	LLWKRAHAESPDCVDNLKSD IGDKGKNKDEGDVCNHEKKT ADLAPHPEEKKKKRSRGFRDR KVMHEYENRIIRAYSTPDKIFR YFATLKVISEPGAEVFMTP EDFVRSITPNEKQPEHLGLD QYIIKRFDGKKISQEREKFA DEGSIFYTLGE