

Anti-ESM1 antibody (20-184 aa) (STJ11106381)

STJ11106381

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	11082
Gene Symbol	ESM1
Uniprot ID	ESM1_HUMAN
Immunogen	
Immunogen Region	20-184 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 20-184 of human ESM1 (NP_008967.1).
Immunogen Sequence	WSNNYAVDCPQHCDSSSECKS SPRCKRTVLDDCGCCRVCAA GRGETCYRTVSGMDGMKCGP GLRCQPSNGEDPFGEEFGIC KDCPYGTFGMDCRETCNCQS GICDRGTGKCLKPFFQYSV TKSSNRFVSLTEHDMASGDG NIVREEVVKENAAGSPVMRK WLNPR