

Anti-CLASRP antibody (1-300 aa) (STJ11106382)

STJ11106382

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	11129
Gene Symbol	CLASRP
Uniprot ID	CLASR_HUMAN
Immunogen	
Immunogen Region	1-300 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CLASRP (NP_008987.2).
Immunogen Sequence	MWHEARKHERKLRGMMVDYK KRAERRREYYEKIKKDPQF LQVHGRACKVHLDSAVALAA ESPVNMMPWQGDNNMIDRF DVRALHDHIPDYTPLLTTI SPEQESDERKONYERYGLV QNDFAGISEEQCLYQIYIDE LYGGLQRPSEDEKKKLAEEK ASIGYTYEDSTVAEVEKAAE KPEEEESAAEEESNSDEDEV IPDIDVEVDVDELNQE QVAD LNKQATTYGMADGDFVRML