

Anti-ANKRD6 antibody (423-722 aa) (STJ11106384)

STJ11106384

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	22881
Gene Symbol	ANKRD6
Uniprot ID	ANKR6_HUMAN
Immunogen	
Immunogen Region	423-722 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 423-722 of human ANKRD6 (NP_055757.3).
Immunogen Sequence	NQLEATVEEIKAEIGSVQDK MNTKLGQMENKTQHQMRLVD KLMVERLSAERTECLNRQQ HSDTEKHGEKQISLVDEL KTWCMLKIQNLEQKLSGDSR ACRAKSTPSTCVDQLVVTAG PAAASDSSPPVVRPKEKALN STATQRLQELSSSDCTGSR LRNVKVQTALLPMNEAARSD QQAGPCVNRGTQTKKSGKSG PTRHRAQQPAASSTCGQPPP ATGSEQTGPHIRDTSQL