

Anti-PIP5K1C antibody (450-550 aa) (STJ11106390)

STJ11106390

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	23396
Gene Symbol	PIP5K1C
Uniprot ID	PI51C_HUMAN
Immunogen	
Immunogen	450-550 aa
Region	
Specificity	A synthetic peptide corresponding to a sequence within amino acids 450-550 of human PIP5K1C (NP_036530.1).
Immunogen	LKSSPSKKGRGGALLAVKPL GPTAAFSASQIPSEREEAQY DLRGARSYPTLEDEGRPDLL PCTPPSFEEATTASIATTLS
Sequence	STSLSIPERSPSETSEQPRY R