

Anti-CCDC9 antibody (100-200 aa) (STJ11106392)

STJ11106392

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	26093
Gene Symbol	CCDC9
Uniprot ID	CCDC9_HUMAN
Immunogen	
Immunogen Region	100-200 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human CCDC9 (NP_056418.1).
Immunogen Sequence	GGRAGMGRASRSWEGSPGEQ PRGGGAGGRGRGRGRGSPH LSGAGDTSISDRKSKEWEER RRQNIKMNEMEKEIAEYER NQREGVLEPNPVRNFLDDPR R