

Anti-CTDSPL2 antibody (325-466 aa) (STJ11106394)

STJ11106394

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	51496
Gene Symbol	CTDSPL2
Uniprot ID	CTSL2_HUMAN
Immunogen	
Immunogen Region	325-466 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 325-466 of human CTDSPL2 (NP_057480.2).
Immunogen Sequence	YVRLRPFFREFLERMSQMYE IILFTASKKVYADKLLNILD PKKQLVRHRLFREHCVCVQG NYIKDLNILGRDLSKTIID
Sequence	NSPQAFAYQLSNGPIESWF MDKNDNELLKLIPFLEKIVE LNEDVRPHIRDRLHDLPL PD