

Anti-TRMT1L antibody (1-300 aa) (STJ11106402)

STJ11106402

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	81627
Gene Symbol	TRMT1L
Uniprot ID	TRMT1L_HUMAN
Immunogen	
Immunogen Region	1-300 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TRMT1L (NP_112196.3).
Immunogen Sequence	MENMAEEELLPLEKEEVEVA QVQVTPARDSAGVPAPAPD SALDSAPTPASAPAPAPALA QAPALSPSLASAPEEAKSKR HISIQRLADLENLAFVTDG NFDSASSLNSDNLDAQNRQA CPLCPKEKFRACNSHKLRRLH LQNLHWKVSVEFEGYRMCIC HLPCRVPKPNIIGEQITSKM GAHYHCIICSATITRTDML GHVRRHMNKGETKSSYIAAS TAKPPKEILKEADTDVQVC