

Anti-TRIM7 antibody (212-511 aa) (STJ11106403)

STJ11106403

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	81786
Gene Symbol	TRIM7
Uniprot ID	TRIM7_HUMAN
Immunogen	
Immunogen Region	212-511 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 212-511 of human TRIM7 (NP_976038.1).
Immunogen Sequence	EQEKVGAEFQALRAFLVEQE GRLLGRLEELSREVAQKQNE NLAQLGVEITQLSKLSSQIQ ETAQKPDLDLQEFKSTLSR CSNVPGPKPTTVSSEMKNKV WNVSLKTFVLKGMLKKFKED LRGELEKEEKVELTDPDTA NPRILSLDLKGVRGGERAQ DLPNHPCRFDNTNTRVLASCG FSSGRHHWEVEVGSKDGWAF GVARESRRKGLTPFTPPEG VWALQLNGGQYWAVTSPER