

Anti-SLC25A2 antibody (21-101 aa) (STJ11106405)

STJ11106405

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	83884
Gene Symbol	SLC25A2
Uniprot ID	ORNT2_HUMAN
Immunogen	
Immunogen Region	21-101 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 21-101 of human SLC25A2 (NP_114153.1).
Immunogen Sequence	TACVLTGQPFDTIKVKMQTF PDLYKGLTDCFLKTYAQVGL RGFYKGTGPALMAYVAENSV LFMCYGFCQQFVRKVAGMDK Q