

Anti-PAXBP1 antibody (200-300 aa) (STJ11106408)

STJ11106408

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	94104
Gene Symbol	PAXBP1
Uniprot ID	PAXB1_HUMAN
Immunogen	
Immunogen Region	200-300 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human PAXBP1 (NP_057715.2).
Immunogen Sequence	TGGAFSNALSSLNVLRPGEI PDAAFIHAARKKKQMARELG DFTPHDNEPGKGRLVREDEN DASDDEDDDEKRRIVFSVKE KSQRQKIAEEIGIEGSDDDA L