

Anti-TBC1D16 antibody (1-200 aa) (STJ11106409)
STJ11106409

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	125058
Gene Symbol	TBC1D16
Uniprot ID	TBC16_HUMAN
Immunogen	
Immunogen Region	1-200 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human TBC1D16 (NP_061893.2).
Immunogen Sequence	MSLGRLLRRASSKASDLLTL TPGGSGSGSPSVLDGEIYS KNNVCVHPPEGLQGLGEHHP GYLCLYMEKDEMLGATLILA WVPNSRIQRQDEEALRYITP ESSPVRKAPRPRGRTRSSG ASHQPSPTELRPTLTPKDED ILWAQSVDPDRMLASPAPED EEKLAQGLGVDGAQPASQPA CSPSGILSTVSPQDVTEEGR