

Anti-ADAMTS9 antibody (1520-1750 aa) (STJ11106425)
STJ11106425

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	56999
Gene Symbol	ADAMTS9
Uniprot ID	ATS9_HUMAN
Immunogen	
Immunogen Region	1520-1750 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1520-1750 of human ADAMTS9 (NP_891550.1).
Immunogen Sequence	RHVGCCQIGTHKIARETECNP YTRPESERDCQGPRCPLYTW RAEWQECTKTCGEGSRYRK VVCVDDNKNEVHGARCDVSK RPVDRECSLQPC EYVWITG EWSECSVTCGKG YKQRLVSC SEIYTGKENYEYSYQTTINC PGTQPPSVHPCYL RDCPVSA TWRVGNWGS CSVSCGVGMQ RSVQCLTNEDQPSHLCHTDL KPEERKTCRNVYNCELPQNC KEVKRLKGASE