

Anti-Ubiquitinyl-Histone H2B-K120 antibody (STJ11106430)

STJ11106430

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat/Other

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	3017/8339/8343/8344/8346/8347 8349
Gene Symbol	H2BC4.H2BC6.H2BC7.H2BC8.H2BC10 H2BC21
Uniprot ID	H2B1C_HUMAN H2B2E_HUMAN
Immunogen	VTKYT
Immunogen Region	
Specificity	A synthetic ubiquitinated peptide around K120 of human HIST1H2BC (NP_066407.1).
Immunogen Sequence	VTKYT