

Anti-KRT23 antibody (1-422 aa) (STJ11106435)

STJ11106435

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	25984
Gene Symbol	KRT23
Uniprot ID	K1C23_HUMAN
Immunogen	
Immunogen Region	1-422 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-422 of human KRT23 (NP_056330.3).
Immunogen Sequence	MNSGHSFSQTPSASFHGAGG GWGRPRSFPRAPTVHGGAGG ARISLSFTTRSCPPPGGSWG SGRSSPLLGGNGKATMQNLN DRLASYLEKVRALLEEANKML ESRILKWHQQRDPGSKKDYS QYEENITHLQEQIVDGKMTN AQIILLIDNARMAVDDFNLK YNEHSFKKDLIEIEGLRR TLDNLTIVTTDLQEQVEGMR KELILMKKHHEQEMEKKHHVP SDFNVNVKVDGTPREDLIK