

Anti-Cre antibody (1-343 aa) (STJ11106443)
STJ11106443

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	All

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen	1-343 aa
Region	
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-343 of e. coli Cre (WP_097745291.1)
Immunogen	MSNLLTVHQNLPALPIDATS DEVRKNLMDMFRDRQAFSEH TWKMLLSVCRSWAAWCKLNN RKWFPAEPEDVRDYLLYLQA
Sequence	RGLAVKTIQQHLGQLNMLHR RSGLP RPSPDSNAVSLVMRRI RKENVDAGERAKQALAFERT DFDQVRSLMENS DRCQDIRN LAFLGIAYNTLLRIAEIARI RVKDISRTDGGRM LIHIGRT KTLVSTAGVEKALSLGVTKL VERWISVSGVADDPNNYLF