

Anti-Tulp2 antibody (78-305 aa) (STJ11106455)
STJ11106455

GENERAL INFORMATION




Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	56734
Gene Symbol	Tulp2
Uniprot ID	TULP2_MOUSE
Immunogen	
Immunogen Region	78-305 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 78-305 of mouse Tulp2 (NP_032833.2).
Immunogen Sequence	FQSDSSWGLGVGSPFLQENV PQAHLPSGAHSALVTMSYVA DSGSERAPLLSPRGAVYTRG NGPAVRHHLCLWLPDSSSDSDV EEVTMEDIPVISRPPQTNLA NLRGWLASPGPGISQEEKE EEVGSTDARVEDKTPSPDPD PDPTVNSDGDHGLAPCKVE ENTAQKNTETASGIGDEDRE KGEVTESTETNYAPVASKVL QGDDGDASNHNWMTCPQP RIPGPRLG



Immunohistochemistry analysis of paraffin-embedded Rat testis using Tulp2 Rabbit pAb (STJ11106455) at dilution of 1:200 (40x lens). High pressure antigen retrieval performed with 0.01M Citrate Buffer (pH 6.0) prior to immunohistochemistry staining.



Immunohistochemistry analysis of paraffin-embedded Mouse testis using Tulp2 Rabbit pAb (STJ11106455) at dilution of 1:200 (40x lens). High pressure antigen retrieval performed with 0.01M Citrate Buffer (pH 6.0) prior to immunohistochemistry staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081