

Anti-TAPBPL antibody (19-170 aa) (STJ11106464)

STJ11106464

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	55080
Gene Symbol	TAPBPL
Uniprot ID	TPSNR_HUMAN
Immunogen	
Immunogen Region	19-170 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 19-170 of human TAPBPL (NP_060479.3).
Immunogen Sequence	AETKPHPAEGQWRAVDVLD CFLAKDGAHRGALASSEDRA RASLVLKQVPVLDDGSLEDF TDFQGGTLAQDDPPIIFEAS VDLVQIPQAEALLHADCSGK EVTCEISRYFLQMTETTVKT AAWFMANMQVSGGGPSISLV MKTPRVTKNEAL