

Anti-PCDH7 antibody (28-300 aa) (STJ11106470)

STJ11106470

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	5099
Gene Symbol	PCDH7
Uniprot ID	PCDH7_HUMAN
Immunogen	
Immunogen Region	28-300 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 28-300 of human PCDH7 (NP_001166994.1).
Immunogen Sequence	AAAKQLRLRYRLAEEGPADVRIGNVASDLGIVTGSGEVTFSL ESGSEYLKIDNLGTGELSTSERRIDREKL PCCQMIFDENE CFLDFEVSIVIGSQSWVDLFEGQVIVLDINDNPTTFPSPVLTLTVEENRPVGTLYLLPTATDRDFGRNGIER YELLQEPG GGSGGGESRRAGAADSAPYPGGGGNGASGGGSGGSKRRILDASEGGGGTNPGGRSSVFELQVADTPDGEKQPQLIVKGAL