

Anti-RECQL5 antibody (400-500 aa) (STJ11106473)

STJ11106473

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	9400
Gene Symbol	RECQL5
Uniprot ID	RECQ5_HUMAN
Immunogen	
Immunogen Region	400-500 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 400-500 of human RECQL5 (NP_004250.4).
Immunogen Sequence	DALVTFCEELGCRHAAIAKY FGDALPACAKGCDHCQNPTA VRRRLEALERSSSSWSKTCIG PSQGNFGDPPELYEGGRKGYG DFSRYDEGSGGSGDEGRDEA H