

Anti-GTF3C1 antibody (1600-1700 aa) (STJ11106476)

STJ11106476

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	
Instruction	

TARGET INFORMATION

Gene ID	2975
Gene Symbol	GTF3C1
Uniprot ID	TF3C1_HUMAN
Immunogen	
Immunogen Region	1600-1700 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1600-1700 of human GTF3C1 (NP_001273171.1).
Immunogen Sequence	NEVIKSLGKDGSLDEDEEE DDLDEGVGGKRRSMEVKPAQ ASHTNYLLMRGYSPGIVST RNLNPNDIVVNSCQMKFQL RCTPVPARLRPAAAPLEELT M