

Anti-PIGF antibody (1-100 aa) (STJ11106479)

STJ11106479

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	5281
Gene Symbol	PIGF
Uniprot ID	PIGF_HUMAN
Immunogen	
Immunogen Region	1-100 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PIGF (NP_002634.1).
Immunogen Sequence	MKDNDIKRLLYTHLLCIFSI ILSVFIPSLFLENFSILETH LTWLCICSGFVTAVNLVLYL VVKPNTSSKRSSLSHKVTGF LKCCYFLMSCFSFHVIFVL

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081