

Anti-PSKH1 antibody (295-424 aa) (STJ11106480)

STJ11106480

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	5681
Gene Symbol	PSKH1
Uniprot ID	KPSH1_HUMAN
Immunogen	
Immunogen Region	295-424 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 295-424 of human PSKH1 (NP_006733.1).
Immunogen Sequence	MPFEDDNRTLYRQILRGKY SYSGEPWPSVSNLAKDFIDR LLTVDPGARM TALQALRHPW VVSM AASSSMKNLHRSISQN LLKRASSRCQSTKSAQSTRS SRSTRSNKSRVRERELREL NLR YQQQYNG