

Anti-SNX2 antibody (1-130 aa) (STJ11106482)

STJ11106482

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC-P: 1:50-1:200 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	

TARGET INFORMATION

Gene ID	6643
Gene Symbol	SNX2
Uniprot ID	SNX2_HUMAN
Immunogen	
Immunogen Region	1-130 aa
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human SNX2 (NP_003091.2).
Immunogen Sequence	MAAEREPPPLGDGKPTDFED LEDGEDLFTSTVSTLESSPS SPEPASLPAEDISANSNGPK PTEVVLDDDDREDLFAEATEE VSLDSPEREPILSSEPSPAV TPVTPPTTLAPRIESKSMSA PVIFDRSREE